

BIOCLEAN CHEMPREP

Sterile Prep Mat S-BCPM

About

The BioClean ChemPrep™ Prep Mat enhances workbench safety for technicians, nurses and pharmacists when handling or administering chemotherapy drugs. The Prep Mat is a highly-absorbent, three-layered preparation mat consisting of two layers of polyethylene with a central layer of air-laid paper. The upper layer allows even distribution of liquid, the middle layer gives absorption properties and the bottom layer provides a defence against leakages.

Specifications

MATERIAL:

Top layer - Polyolefin apertured film

Middle layer - Air-laid paper, Bottom

Layer - Polyethylene

WEIGHT:

210gsm

LENGTH:

56cm (22") or 40cm (16")

COMPATIBILITY:

ISO Class 5 & EU GMP Grade A

COLOUR:

Blue

STERILIZATION:

Gamma irradiation

Features

- Latex-free 3-layer construction
- Mesh upper layer for even distribution
- Highly absorbent middle layer
- Impermeable bottom layer

HIGHLY ABSORBENT

LOW-LINTING

3-LAYER CONSTRUCTION



Quality Standards

- Manufactured in a facility holding ISO 9001:2015
- Processed in an NEBB certified ISO Class 5 environment

Sterilization

Method: Gamma irradiation
Minimum Dose: 25kGy
Sterility Assurance Level: 10^{-6}

Shelf Life & Storage

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

Physcial/General Properties

Property Measured	Method Used	Units of Measurement	Typical Value
Sorptive Capacity	IEST-RP-CC004.3	ml/m ²	852
Sorptive Capacity	IEST-RP-CC004.3	ml/prep mat	83.3



TO ORDER

Re-order code

S-BCPM-4025B
S-BCPM-5640B

Size

40cm x 25cm
56cm x 40cm

Packing

S-BPCM-4025B:

One piece per sealed inner PE bag; one inner bag per second PE bag; 50 double bags per sealed outer PE bag; four outer bags per lined carton (200 pieces)

S-BPCM-5640B:

One piece per sealed inner PE bag; one inner bag per second PE bag; 25 double bags per sealed outer PE bag; four outer bags per lined carton (100 pieces)

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd

Ground Floor, 15 Kings
Court
Willie Snaith Road,
Newmarket,
Suffolk, CB8 7SG
UK
Tel: +44 1638 663338
S-BCPM/1811/PDS20

Ansell US

111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2020 Ansell Limited. All Rights Reserved.