

Excell BEXS

Sterile Disposable Nitrile Cleanroom Glove



Accelerator-free sterile nitrile gloves, offering users enhanced protection

- Reduced allergy risks: BioClean[™] Excell BEXS sterile nitrile gloves are powderfree, latex-free and chemical accelerator-free, eliminating Type I latex and Type IV chemical allergies and sensitivities
- Clean processed: They undergo deionized washing to reduce surface particulates, ionic extractables and nonvolatile residues
- More sustainable packaging: Packed in recyclable* plastic packaging and delivered in recycled cardboard shipper cases
- Decreased contamination risks: Nonparticulating packaging reduces potential contaminant hazards in the cleanroom environment
 - *Always check the recyclable status of HDPE locally, as this plastic material may not be considered suitable for recycling in your country
- Chemical resistance: Protection from chemicals and chemotherapy drugs
- Hand-specific ergonomic design: Improved dexterity for skilled handling
- Gamma irradiation sterilized: Sterility assurance level (SAL) of 10-6

Industries

- Sterile Processing Department
- Biotechnology Manufacturing
- Controlled and Critical Environments
- Laboratory Research and Development
- Medical Device Manufacturing
- Pharmaceutical Manufacturing







Excell BEXS

Sterile Disposable Nitrile Cleanroom
Glove

TECHNICAL DATA SHEET

| | Product Information |
|--------------------------------------|--|
| Material | Nitrile |
| Color | White |
| Shape | Hand Specific |
| Cuff | Beaded |
| Manufacturing/QMS Audit Standards | Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D |
| Regulatory | CE 0598, EN 420:2003 + A1:2009, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Food Contact, Category III |
| Packaging | 1 pair per inner wrap; 1 inner wrap per pouch; 10 pouches per outer bag; 20 outer bags per lined carton (200 pairs) |
| Storage | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Country of Origin | Malaysia |
| Available sizes | 6.0, 6.5, 7 0, 7.5, 8.0, 8.5, 9.0 |
| External Glove Surface | Textured Fingers and Palm |
| Internal Glove Surface | Chlorinated |
| Sterilization Method | GAMMA irradiation (25 kGy) |
| Sterilization Minimum Dose | 25kGy |
| Sterility Assurance Level | 10-6 |
| Cleanroom Class | Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms |
| Shelf Life | Five (5) years from date of manufacture |
| Protein Level | N/A: contains no natural rubber latex |
| Anti-static | Yes |





Excell BEXS

Sterile Disposable Nitrile Cleanroom Glove

| Physical Properties | | | | | | | | Testing Method | |
|--|---------------------|----------|--------|-----|-----|----------|-----------------|----------------|--|
| Sizes | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 | | |
| Typical Length (mm/in) | | 300 / 12 | | | | | | EN 420 | |
| Palm Width (mm/in) | 79 | 84 | 89/3.5 | 95 | 101 | 108/4.2 | 115 | | |
| Freedom from Holes | 1.5 AQL Performance | | | | | | EN 374-2 | | |
| Typical Particle Count ≥0.5µm (counts / cm²) | <1200 | | | | | | IEST-RP-CC005.4 | | |
| Target Single Wall Palm Thickness (mm/mil) | 0.14 / 5.51 | | | | | | EN 455-2 | | |
| Target Single Wall Finger Thickness (mm/mil) | 0.17 / 6.69 | | | | | | EN 455-2 | | |
| Target Single Wall Cuff Thickness (mm/mil) | 0.11 / 4.33 | | | | | EN 455-2 | | | |
| Ultimate tensile strength (MPa) Before Aging | Min. 15 | | | | | | ASTM D412-06a | | |

IONIC CONTENT

| Concentration in µg/cm² | Typical | Concentration in µg/cm² | Typical | |
|-------------------------|-------------------|-------------------------|-------------------|--|
| Ammonium | Available on File | Nitrate | <0.25 | |
| Bromide | <0.01 | Nitrite | Available on File | |
| Calcium | <0.55 | Phosphate | <0.01 | |
| Chloride | <0.35 | Potassium | <0.10 | |
| Fluoride | < 0.01 | Sodium | < 0.05 | |
| Lithium | < 0.01 | Sulphate | <0.05 | |
| Magnesium | < 0.01 | Zinc | < 0.05 | |

ORDERING INFORMATION

| SIZE | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 |
|----------------|---------|---------|---------|---------|---------|---------|---------|
| REORDER NO. | BEXS-60 | BEXS-65 | BEXS-70 | BEXS-75 | BEXS-80 | BEXS-85 | BEXS-90 |

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE











For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV T: +32 (0) 2 528 74 00 F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center T: +603 8310 6688 F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC US T: +1800 800 0444 US F: +1800 800 0445 CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V. T: +52 442 248 1544 / 248 3133

Ansell Limited T: +61 1800 337 041 F: +61 1800 803 578

UK

Ansell Nitritex T: +44 1638 663338 F: +44 1638 668890

Ansell, $^{\textcircled{R}}$ and $^{\textcircled{M}}$ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking $^{\textcircled{M}}$ 2022 Ansell Limited. All Rights Reserved.

