### DATASHEET



## **Features**

- » High sorption capacity
- » Unique fabric structure
- » Gamma irradiated in accordance with ANSI/AAMI/ISO 11137
- » Certificates of Processing and Compliance are attached to each case
- » Quarterly microbial testing by Agarose Overlay and Limulus Amebocyte Lysate (LAL) endotoxin testing

# **Vertex**<sup>®</sup> Sterile High-Sorption Wipers

TX3042 and TX3049 are 100% continuous-filament, polyester knit fabric wipers made in the USA from virgin stock and processed on Texwipe's revolutionary Vertex<sup>®</sup> manufacturing system. The wiper's unique design provides a very soft feel and high sorptive capacity.

The Vertex<sup>®</sup> manufacturing process cleans and converts wipers in an integrated, hands-free micro-environment. Vertex<sup>®</sup> wipers are manufactured in an ISO Class 3 environment with a small physical footprint which provides for a more sustainable process.

The Vertex<sup>®</sup> system incorporates real time process monitoring and alarms for critical parameters to produce the most consistent cleanroom wipers available.

Quarterly Agar Overlay and LAL testing assures that wipers do not contribute cytotoxic leachables and are nonprygenic.

#### **Benefits**

- » Reduces number of wipers needed for spill control
- » Soft feel for use with delicate surfaces or applications
- » Meets Sterility Assurance Level of 10<sup>-6</sup>
- » Confirms radiation dose met specified dose at time of receiving
- » Meets FDA Guideline for medical device bacterial endotoxin limit of not more than 20 EU/device

#### **Products**

Number	Dry Wiper	s	Packaging					
TX3042	Vertex®	12" x 12" nominal (31 cm x 31 cm) sealed-edge wipers, triple-bagged	100 wipers/bag (5 inner bags of 20 wipers) 5 bags/case 100 wipers/bag (5 inner bags of 20 wipers) 5 bags/case					
TX3049	Vertex®	9" x 9" nominal (23 cm x 23 cm) sealed-edge wipers, triple-bagged						
Number	Pre-wett	ed Wipers	Packaging					
TX3042P	Vertex®	12" x 12" nominal (31 cm x 31 cm) sealed-edge wipers pre-wetted with 70% USP-Grade IPA / 30% USP Purified Water, triple-bagged	25 wipers/bag, 5 bags/case					
TX3049P	Vertex®	9" x 9" nominal (23 cm x 23 cm) sealed-edge wipers pre-wetted with 70% USP-Grade IPA / 30% USP Purified Water, triple-bagged	25 wipers/bag, 5 bags/case					

ITW Texwipe is the only wiper company to be ISO 9001:2000, 13485:2003,14001:2004 and OHSAS 18001:1999 registered.



**EXVIPE**<sup>®</sup> For more information contact us at the numbers below or www.texwipe.com

©2012 ITW Texwipe Printed in USA • Effective: May 2012

North America Office P: (800) 839-9473 F: (336) 996-2297 Asia/Pacific Office P: +65 6468 9433 F: +65 6468 6772 Europe/Middle East Office P: +31 88 1307 400 F: +31 88 1307 499

# **Vertex® Dry Wipers**

## **Performance and Contamination Characteristics**

Property	<b>Typical Value</b>	Test Method*	Property	Typical Value	Test Method*
Absorbency Sorptive capacity	450 mL/m <sup>2</sup>	1,TM20	LPC ≥0.5 μm	9.5 x 10 <sup>6</sup> particles/m <sup>2</sup>	1,TM22
Sorptive rate	0.3 seconds	1,TM20	Particles and fibers Particles 0.5–5.0 µm 5.0–100 µm Fibers >100 µm	3.0 x 10 <sup>6</sup> particles/m <sup>2</sup> 66,000 particles/m <sup>2</sup> 200 fibers/m <sup>2</sup>	1,2,TM22 1,2,TM22 1,2,TM22
			Nonvolatile residue IPA extractant DIW extractant	0.02 g/m <sup>2</sup> 0.01 g/m <sup>2</sup>	1,TM1 1,TM1

**Note:** The data in these tables represent typical analyses of these wipers at the time of publication. These are not specifications.

#### \*Test Methods

2

1 - IEST-RP-CC004.3, 2004; www.iest.org.

- ASTM E2090-06, 2006; E2090-06, ASTM International; www.astm.org.
- TM Texwipe Test Method available upon request. Visit www.texwipe.com or contact Texwipe Customer Service at + 1-336-996-7046.

# **Wiper Consistency Data**

