

3M Products for Environmental Sample Collection


Collect ▶

Transport ▶

Prepare ▶

Analyze ▶

Description	Catalog Number	Case Quantity
3M™ Sample Handling		
3M™ Swab-Sampler w/10mL Buffered Peptone Water	RS96010BPW	100
3M™ Swab-Sampler w/10mL D/E Neutralizing Broth	RS96010DE	100
3M™ Swab-Sampler w/10mL Lethen Broth	RS96010LET	100
3M™ Swab-Sampler w/10mL Neutralizing Buffer	RS96010NB	100
3M™ Swab-Sampler w/1mL Lethen Broth	RS9601LET	100
3M™ Swab-Sampler w/4mL Lethen Broth	RS9604LET	100
3M™ Swab-Sampler w/4mL Neutralizing Buffer	RS9604NB	100
3M™ Swab-Sampler w/5mL Lethen Broth	RS9605LET	100
3M™ Enviro Swab	ENVSWB25	25
3M™ Enviro Swab, Dry	ENVSWB25D	25
3M™ Sponge-Stick, 1.5" x 3", dry, 30 oz. bag	SSL100	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 30 oz. bag, 5 bags of 20	SSL10BPW	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Lethen Broth, 5 bags of 20	SSL10LET	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Lethen Broth, 30 oz. bag, 2 gloves, 5 bags of 20	SSL10LETG	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 30 oz. bag, 5 bags of 20	SSL10NB	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 30 oz. bag, 2 gloves, 5 bags of 20	SSL10NBG	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 6" x 9" bag, 5 bags of 20	SSL10NB6X9	100
3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL D/E Broth, 30 oz. bag, 5 bags of 20	SSL10DE	100
3M™ Hydrated Sponge, 1.5" x 3" with Buffered Peptone Water Broth, 5 bags of 20	HS10BPW	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 2 gloves, 5 bags of 20	HS10BPW2G	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Lethen Broth, 5 bags of 20	HS10LET	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Lethen Broth, 2 gloves, 5 bags of 20	HS10LET2G	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 5 bags of 20	HS10NB	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 2 gloves, 5 bags of 20	HS10NB2G	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 50 oz. bag	HS119510NB	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 24 oz. bag, 2 gloves	HS2410NB2G	100
3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL D/E Neutralizing Buffer, 2 gloves	HS10DE2G	100
3M™ Clean-Trace Hygiene Monitoring Systems		
3M™ Clean-Trace™ Surface ATP Swab	UXL100	100
3M™ Clean-Trace™ Surface Protein Plus Swab	PR0100	100
3M™ Clean-Trace™ Surface Protein (Allergen) Swab	ALLTEC60	60



Food Safety

3M Center, Building 275-5W-05
St. Paul, MN 55144-1000
U.S.A.
1-800-328-6553
www.3M.com/foodsafety

Please recycle. Printed in U.S.A.
© 3M 2012. All Rights Reserved.
70-2011-5021-9

3M and Clean-Trace are trademarks of 3M Company.
Used under license in Canada.